

Product datasheet for **HP200953**

NKIRAS2 Human qPCR Primer Pair (NM_001001349)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	28511
Forward Sequence:	GGACAAGAAGGAGGTCACCATC
Reverse Sequence:	GCTTCACCTTCTCTGACTTGGC
ACCN:	<u>NM_001001349</u> , <u>NM_001001349.1</u> , <u>NM_001001349.2</u> , <u>BC063498</u> , <u>BC063498.1</u> , <u>BC007316</u> , <u>BC007450</u> , <u>BC009059</u> , <u>BG717103</u>
UniProt ID:	<u>Q9NYR9</u>
Synonyms:	kappaB-Ras2; KBRAS2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)