

Product datasheet for **HP200955**

LIM domain only 3 (LMO3) Human qPCR Primer Pair (NM_001001395)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55885
Forward Sequence:	ATGGTGATGCGTGCCAAGGACA
Reverse Sequence:	AACCTTCCTCGTAGTCCGCTG
ACCN:	NM_001001395 , NM_001001395.1 , NM_001001395.2 , BC026311 , BC026311.1 , BC050085 , BC017777 , BG714670 , BI553855 , BQ029059
UniProt ID:	Q8TAP4
Synonyms:	RBTN3; RBTNL2; Rhom-3; RHOM3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)