

Product datasheet for **HP200979**

HRP2 (HDGFRP2) Human qPCR Primer Pair (NM_001001520)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84717
Forward Sequence:	CAAGAGTAGCGACAACAGTGGC
Reverse Sequence:	GCTTTCCGAGTTCTCCTCTCG
ACCN:	<u>NM_001001520</u> , <u>NM_001001520.1</u> , <u>NM_001001520.2</u> , <u>BC000755</u> , <u>BC009449</u> , <u>BE745571</u> , <u>BI826321</u> , <u>BI905780</u> , <u>BM147540</u> , <u>BU738961</u> , <u>NM_001001520.3</u>
UniProt ID:	<u>Q7Z4V5</u>
Synonyms:	HDGF-2; HDGF2; HDGFRP2; HRP-2; HRP2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)