

## Product datasheet for **HP201048**

### LCN10 Human qPCR Primer Pair (NM\_001001712)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	414332
Forward Sequence:	GTTTGGGAACGCCTGTGCATAC
Reverse Sequence:	GGTAGACCAAGCCGTAGCTGTA
ACCN:	<u><a href="#">NM_001001712</a></u> , <u><a href="#">NM_001001712.1</a></u> , <u><a href="#">NM_001001712.2</a></u> , <u><a href="#">BC133045</a></u> , <u><a href="#">BC036875</a></u> , <u><a href="#">BC040963</a></u> , <u><a href="#">BC046377</a></u> , <u><a href="#">BC064596</a></u> , <u><a href="#">BC144440</a></u> , <u><a href="#">BC144441</a></u> , <u><a href="#">BC171741</a></u> , <u><a href="#">NM_001001712.3</a></u>
UniProt ID:	<u><a href="#">Q6JVE6</a></u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)