

Product datasheet for **HP201163**

GATA3 Human qPCR Primer Pair (NM_001002295)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2625
Forward Sequence:	ACCACAACCACACTCTGGAGGA
Reverse Sequence:	TCGGTTTCTGGTCTGGATGCCT
ACCN:	BC006793 , NM_001002295 , NM_001002295.1 , BC006793.1 , BC003070 , BC006839 , BG287052 , BM152193 , NM_001002295.2
UniProt ID:	P23771
Synonyms:	HDR; HDRS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)