

Product datasheet for **HP201167**

C10orf78 (SFR1) Human qPCR Primer Pair (NM_001002759)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	119392
Forward Sequence:	CCTCTGCGAATCCATCATCTCC
Reverse Sequence:	GACGTTTCACCACATTGTAAGAGG
ACCN:	<u>NM_001002759</u> , <u>NM_001002759.1</u> , <u>BC020892</u> , <u>BC043256</u> , <u>BC047587</u> , <u>BC140835</u> , <u>BE960638</u> , <u>BG942691</u>
UniProt ID:	<u>Q86XK3</u>
Synonyms:	bA373N18.1; C10orf78; MEI5; MEIR5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)