

Product datasheet for **HP201180**

CHTF8 Human qPCR Primer Pair (NM_001002847)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54921
Forward Sequence:	ACCAGGTCCTATGTCCAACCCA
Reverse Sequence:	GCTCTCAGGTTAGGTCCAGATC
ACCN:	BC018700 , BC018700.2 , BC069059 , BM820172 , BP195395 , BP222147 , BU631276 , BU739855
UniProt ID:	P0CG13
Synonyms:	CHTF8
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)