

Product datasheet for **HP201194**

GPATCH8 Human qPCR Primer Pair (NM_001002909)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23131
Forward Sequence:	GCTGAAGATGCTACCGAACGGC
Reverse Sequence:	CTGAGGTCTTCCAAGGCTTTGG
ACCN:	<u>NM_001002909</u> , <u>NM_001002909.1</u> , <u>NM_001002909.2</u> , <u>NM_001002909.3</u> , <u>BC017808</u> , <u>BC024147</u> , <u>BC063601</u> , <u>BC136629</u> , <u>BC140010</u> , <u>BC144323</u> , <u>BM722954</u> , <u>NM_001002909.4</u>
UniProt ID:	<u>Q9UKJ3</u>
Synonyms:	GPATC8; KIAA0553
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)