

Product datasheet for **HP201197**

C1orf173 (ERICH3) Human qPCR Primer Pair (NM_001002912)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	127254
Forward Sequence:	GAAAGCACTGCTCAGGTGAAGG
Reverse Sequence:	GCACCAAGAGCCAGGACATAGG
ACCN:	<u>NM_001002912</u> , <u>NM_001002912.1</u> , <u>NM_001002912.2</u> , <u>NM_001002912.3</u> , <u>NM_001002912.4</u> , <u>BC035114</u> , <u>BC044605</u> , <u>BC073916</u> , <u>BI832067</u>
UniProt ID:	<u>Q5RHP9</u>
Synonyms:	C1orf173
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)