

Product datasheet for **HP201210**

POLR1F Human qPCR Primer Pair (NM_001002926)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	221830
Forward Sequence:	CTGCTGGAGTATTCTGCATTCGG
Reverse Sequence:	GCAGCTTCCTCAGTGCCATTTTC
ACCN:	<u>NM_001002926</u> , <u>NM_001002926.1</u> , <u>BC130298</u> , <u>BC013179</u> , <u>BC014574</u> , <u>BC056267</u> , <u>BC070358</u> , <u>BC130300</u> , <u>NM_001002926.2</u>
UniProt ID:	<u>Q3B726</u>
Synonyms:	A43; RPA43; TWISTNB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)