

Product datasheet for **HP201306**

ALG11 Human qPCR Primer Pair (NM_001004127)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	440138
Forward Sequence:	CATCCATACTGCAATGCTGGTGG
Reverse Sequence:	GACCGTTGACATTAACATCGCCG
ACCN:	BC010857 , NM_001004127 , NM_001004127.1 , BC073862 , BC108677 , BC111022 , BC142998 , BC143059 , BC150316 , BP373026
UniProt ID:	Q2TAA5
Synonyms:	CDG1P; GT8
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)