

Product datasheet for **HP201335**

LIN28B Human qPCR Primer Pair (NM_001004317)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	389421
Forward Sequence:	CCTGTTTAGGAAGTGAAAGAAGAC
Reverse Sequence:	CACTTCTTGGCTGAGGAGGTAG
ACCN:	<u>NM_001004317</u> , <u>NM_001004317.1</u> , <u>NM_001004317.2</u> , <u>NM_001004317.3</u> , <u>BC127712</u> , <u>BC127713</u> , <u>BC137526</u> , <u>BC137527</u> , <u>BC141959</u> , <u>BC141960</u> , <u>BC150537</u> , <u>BC150538</u> , <u>NM_001004317.4</u>
UniProt ID:	<u>Q6ZN17</u>
Synonyms:	CSDD2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)