

Product datasheet for **HP201356**

ZYG11A Human qPCR Primer Pair (NM_001004339)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	440590
Forward Sequence:	CTCCACAGGCTGTTACAGAGA
Reverse Sequence:	GCGTGTTAGGTTGAGAGCGCAA
ACCN:	<u>NM_001004339</u> , <u>NM_001004339.1</u> , <u>NM_001004339.2</u> , <u>BC119789</u> , <u>BC119815</u> , <u>NM_001004339.3</u>
UniProt ID:	<u>Q6WRX3</u>
Synonyms:	ZYG11
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)