

## Product datasheet for **HP201567**

### OR4F17 Human qPCR Primer Pair (NM\_001005240)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	81099
Forward Sequence:	TGGGAGTGAGATGGTGATCCTC
Reverse Sequence:	CAACTGGCTCACCGAATGGAGA
ACCN:	<u><a href="#">NM_001005240</a></u> , <u><a href="#">NM_001005240.1</a></u> , <u><a href="#">BC136848</a></u> , <u><a href="#">BC136867</a></u> , <u><a href="#">NM_001005240.2</a></u>
UniProt ID:	<u><a href="#">Q8NGA8</a></u>
Synonyms:	OR4F5; OR4F11P; OR4F18; OR4F19
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)