

Product datasheet for **HP201683**

OR51B5 Human qPCR Primer Pair (NM_001005567)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	282763
Forward Sequence:	AACTCGCCTGTGCTGATACCAC
Reverse Sequence:	TCTGGAGGCAATGCTCAGGACA
ACCN:	<u>NM_001005567</u> , <u>NM_001005567.1</u> , <u>NM_001005567.2</u> , <u>BC022184</u> , <u>BC136707</u> , <u>BC136708</u> , <u>BG493214</u> , <u>BG776801</u>
UniProt ID:	<u>Q9H339</u>
Synonyms:	HOR5'Beta5; OR11-37
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)