

Product datasheet for **HP201815**

SFI1 Human qPCR Primer Pair (NM_001007467)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9814
Forward Sequence:	CTCAGAAGAGGCGGTACAAGCA
Reverse Sequence:	CATCGCCAGAGTCTGTCCATG
ACCN:	<u>NM_001007467</u> , <u>NM_001007467.1</u> , <u>NM_001007467.2</u> , <u>BC021576</u> , <u>BC033613</u> , <u>BC041101</u> , <u>BC046221</u> , <u>BC110814</u> , <u>BC129926</u> , <u>BC129945</u> , <u>NM_001007467.3</u>
UniProt ID:	<u>A8K8P3</u>
Synonyms:	hSfi1p; PISD; PPP1R139
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)