

## Product datasheet for **HP201851**

### TTC6 Human qPCR Primer Pair (NM\_001007795)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	319089
Forward Sequence:	GCCATTGAGGTCTTGGACGGAA
Reverse Sequence:	GAGAAATGGCATCCAAGCCTTC
ACCN:	<a href="#">NM_001007795</a> , <a href="#">BC014342</a> , <a href="#">BC038110</a> , <a href="#">BC103914</a> , <a href="#">BC103915</a> , <a href="#">BC103916</a> , <a href="#">BX161415</a>
Synonyms:	MGC119358; MGC119358, MGC119360, MGC119361; MGC119360; MGC119361; tetratricopeptide repeat domain 6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)