

Product datasheet for **HP201935**

TMCO2 Human qPCR Primer Pair (NM_001008740)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	127391
Forward Sequence:	GCCTTATGGAAACGAAGTATTCAG
Reverse Sequence:	TTCGGAGACCACTTCCTTCTGG
ACCN:	<u>NM_001008740</u> , <u>NM_001008740.1</u> , <u>NM_001008740.2</u> , <u>NM_001008740.3</u> , <u>BC108666</u> , <u>BU567634</u> , <u>BU569147</u> , <u>NM_001008740.4</u>
UniProt ID:	<u>Q7Z6W1</u>
Synonyms:	dj39G22.2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)