

## Product datasheet for **HP202004**

### C12orf29 Human qPCR Primer Pair (NM\_001009894)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	91298
Forward Sequence:	GCCATACCTTTGGGCTCGACTA
Reverse Sequence:	GCTGGTTTGAAGTCCTCCTCAAC
ACCN:	<u><a href="#">NM_001009894</a></u> , <u><a href="#">NM_001009894.1</a></u> , <u><a href="#">NM_001009894.2</a></u> , <u><a href="#">BC092422</a></u> , <u><a href="#">BC092422.1</a></u> , <u><a href="#">BC035136</a></u> , <u><a href="#">BC058014</a></u> , <u><a href="#">BX537868</a></u> , <u><a href="#">BX640968</a></u> , <u><a href="#">BX648026</a></u> , <u><a href="#">NM_001009894.3</a></u>
UniProt ID:	<u><a href="#">Q8N999</a></u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)