

## Product datasheet for **HP202035**

### NRG3 Human qPCR Primer Pair (NM\_001010848)

#### Product data:

|                   |   |
|-------------------|---|
| Product Type:     | qPCR Primer Pairs   |
| Gene ID:          | 10718   |
| Forward Sequence: | TGCCGAGACAAGGACCTTG CAT   |
| Reverse Sequence: | TCACAACGGACTCCTTGGTAGC  |
| ACCN:             | <u>NM_001010848, NM_001010848.1, NM_001010848.2, BC136811, BC143665, NM_001010848.4</u>   |
| UniProt ID:       | <u>P56975</u>   |
| Synonyms:         | HRG3; pro-NRG3  |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.   |
| Stability:        | The primer mix is stable for one year from date of shipping.  |



[View online »](#)