

Product datasheet for **HP202039**

PM20D2 Human qPCR Primer Pair (NM_001010853)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	135293
Forward Sequence:	CACATTGGATCTAATGCCTTGAATC
Reverse Sequence:	ACATCCAGTGCCGTCATTGCCA
ACCN:	<u>NM_001010853</u> , <u>NM_001010853.1</u> , <u>NM_001010853.2</u> , <u>BC035036</u> , <u>BC035036.1</u> , <u>BC045583</u> , <u>BX640865</u> , <u>NM_001010853.3</u>
UniProt ID:	<u>Q8IYS1</u>
Synonyms:	ACY1L2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)