

## Product datasheet for **HP202041**

### PIK3R6 Human qPCR Primer Pair (NM\_001010855)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	146850
Forward Sequence:	CTGGAGGAGATTTACTGCTCGC
Reverse Sequence:	GGTCCACAAGTGGGAGGTGATG
ACCN:	<a href="#">NM_001010855</a> , <a href="#">NM_001010855.1</a> , <a href="#">NM_001010855.2</a> , <a href="#">NM_001010855.3</a> , <a href="#">BC110301</a> , <a href="#">BC125169</a> , <a href="#">BC166645</a> , <a href="#">NM_001010855.4</a>
UniProt ID:	<a href="#">Q5UE93</a>
Synonyms:	C17orf38; HsT41028; p84 PIKAP; p87(PIKAP); p87PIKAP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)