

## Product datasheet for **HP202057**

### TSPO2 Human qPCR Primer Pair (NM\_001010873)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	222642
Forward Sequence:	GCCTCTATGCTGTTTCAGCTCAC
Reverse Sequence:	ATGGGATGCCAGATCAGTGCTG
ACCN:	<u><a href="#">NM_001010873</a></u> , <u><a href="#">NM_001010873.1</a></u> , <u><a href="#">NM_001010873.2</a></u> , <u><a href="#">BC144710</a></u> , <u><a href="#">BC104440</a></u> , <u><a href="#">BC104441</a></u> , <u><a href="#">BC137560</a></u> , <u><a href="#">BC137565</a></u> , <u><a href="#">BU655765</a></u> , <u><a href="#">BU657409</a></u> , <u><a href="#">NM_001010873.3</a></u>
UniProt ID:	<u><a href="#">Q5TGU0</a></u>
Synonyms:	BZRPL1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)