

Product datasheet for **HP202097**

BCL2L15 Human qPCR Primer Pair (NM_001010922)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	440603
Forward Sequence:	ATTGCTGGTCGCCTTCGGATGT
Reverse Sequence:	CCACAGTGCCTGGAGTATAGC
ACCN:	<u>NM_001010922</u> , <u>NM_001010922.1</u> , <u>NM_001010922.2</u> , <u>BC127719</u> , <u>BC127720</u> , <u>BC132904</u> , <u>BC132906</u> , <u>BX116361</u> , <u>NM_001010922.3</u>
UniProt ID:	<u>Q5TBC7</u>
Synonyms:	Bfk; C1orf178
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)