

Product datasheet for **HP202116**

C1orf130 (NCMAP) Human qPCR Primer Pair (NM_001010980)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	400746
Forward Sequence:	TTCACCGTGGTTCTGATCCTGC
Reverse Sequence:	CTGAAGGTCACAGTTGCTGGGT
ACCN:	<u>NM_001010980</u> , <u>NM_001010980.1</u> , <u>NM_001010980.2</u> , <u>NM_001010980.3</u> , <u>NM_001010980.4</u> , <u>BC050474</u> , <u>BC127785</u> , <u>BC147014</u> , <u>BC147015</u> , <u>NM_001010980.5</u>
UniProt ID:	<u>Q5T1S8</u>
Synonyms:	C1orf130; MP11
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)