

Product datasheet for **HP202155**

MAK3 (NAA30) Human qPCR Primer Pair (NM_001011713)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	122830
Forward Sequence:	AAGATCTGTCCGAACCCTACTCC
Reverse Sequence:	GCAAACGATGGCACCTACACAC
ACCN:	<u>NM_001011713</u> , <u>NM_001011713.1</u> , <u>NM_001011713.2</u> , <u>BC048983</u> , <u>BC118589</u> , <u>BC122557</u> , <u>BX648700</u> , <u>NM_001011713.3</u>
UniProt ID:	<u>Q147X3</u>
Synonyms:	C14orf35; MAK3; Mak3p; NAT12; NAT12P
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)