

Product datasheet for **HP202159**

XKR9 Human qPCR Primer Pair (NM_001011720)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	389668
Forward Sequence:	GATAGAGTGACTGATTTGAGCATG
Reverse Sequence:	GCACAGCAAGAGACCATGATGG
ACCN:	<u>NM_001011720</u> , <u>NM_001011720.1</u> , <u>BC111381</u> , <u>BC137104</u> , <u>BC137105</u> , <u>BC151845</u> , <u>NM_001011720.2</u>
UniProt ID:	<u>Q5GH70</u>
Synonyms:	XRG9
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)