

Product datasheet for **HP202175**

C6orf174 (SOGA3) Human qPCR Primer Pair (NM_001012279)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	387104
Forward Sequence:	GCAAGAGGAGATGGAGAAGCTG
Reverse Sequence:	AGTTGACAGGCATCCTCCTCGA
ACCN:	<u>BC020441</u> , <u>NM_001012279</u> , <u>NM_001012279.1</u> , <u>NM_001012279.2</u> , <u>BC127760</u> , <u>BC143043</u> , <u>BC143044</u> , <u>BC156253</u> , <u>BC157034</u> , <u>BE552225</u> , <u>BE779627</u> , <u>NM_001012279.3</u>
UniProt ID:	<u>Q5TF21</u>
Synonyms:	C6orf174; dj403A15.3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)