

## Product datasheet for **HP202210**

### **C9orf117 (CFAP157) Human qPCR Primer Pair (NM\_001012502)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	286207
<b>Forward Sequence:</b>	CAGCAGGAACTGGCTAATGAGC
<b>Reverse Sequence:</b>	ACGTCACTGTCCTCTTCATCGC
<b>ACCN:</b>	<u><a href="#">NM_001012502</a></u> , <u><a href="#">NM_001012502.1</a></u> , <u><a href="#">NM_001012502.2</a></u> , <u><a href="#">BC125172</a></u> , <u><a href="#">BC133027</a></u> , <u><a href="#">BC141809</a></u> , <u><a href="#">NM_001012502.3</a></u>
<b>UniProt ID:</b>	<u><a href="#">Q5JU67</a></u>
<b>Synonyms:</b>	C9orf117
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)