

Product datasheet for **HP202264**

PLAC9 Human qPCR Primer Pair (NM_001012973)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	219348
Forward Sequence:	AGGAGACTCAGCTCAGAGCACA
Reverse Sequence:	CAGGTGATCCACGGTCTTCTCT
ACCN:	NM_001012973 , NM_001012973.1 , BC066348 , BC066348.1 , BC090922 , BC117332 , BC117334 , BU731023 , NM_001012973.3
UniProt ID:	Q5JTB6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)