

Product datasheet for **HP202276**

C1orf31 (COA6) Human qPCR Primer Pair (NM_001012985)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	388753
Forward Sequence:	ATCGCCTGCTTGGTGATTGTGG
Reverse Sequence:	GTTCTCATCTAAACACTTCCAGTAC
ACCN:	<u>NM_001012985</u> , <u>NM_001012985.1</u> , <u>NM_001012985.2</u> , <u>BC116455</u> , <u>BC025793</u> , <u>BC035925</u> , <u>BC052810</u> , <u>BC064905</u> , <u>BG195363</u> , <u>BG204581</u>
UniProt ID:	<u>Q5JTJ3</u>
Synonyms:	C1orf31; CEMCOX4; MC4DN13
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)