

Product datasheet for **HP202320**

NHSL2 Human qPCR Primer Pair (NM_001013627)

Product data:

| | |
|-------------------|---|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 340527 |
| Forward Sequence: | TACTCCAGCCAGTCGGAACAC |
| Reverse Sequence: | GGTAGTGATGACTTCCTCTCAGG |
| ACCN: | <u>NM_001013627</u> , <u>NM_001013627.1</u> , <u>NM_001013627.2</u> , <u>BC033261</u> , <u>BC136756</u> |
| UniProt ID: | <u>Q5HYW2</u> |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec. |
| Storage: | Store at -20 °C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)