

## Product datasheet for **HP202456**

### ODF3B Human qPCR Primer Pair (NM\_001014440)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	440836
Forward Sequence:	CTGGCAGTTTCTTCGAGGACCT
Reverse Sequence:	AAGTCCGCGCCAGAATCGTGAA
ACCN:	<u><a href="#">NM_001014440</a></u> , <u><a href="#">NM_001014440.1</a></u> , <u><a href="#">NM_001014440.2</a></u> , <u><a href="#">NM_001014440.3</a></u> , <u><a href="#">BC127934</a></u> , <u><a href="#">BC127935</a></u> , <u><a href="#">BC172351</a></u> , <u><a href="#">BG149596</a></u> , <u><a href="#">BM973820</a></u> , <u><a href="#">NM_001014440.4</a></u>
UniProt ID:	<u><a href="#">A8MYP8</a></u>
Synonyms:	FAP123; ODF3L3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)