

Product datasheet for **HP202482**

C16orf93 (CCDC189) Human qPCR Primer Pair (NM_001014979)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	90835
Forward Sequence:	ACCGCTACTTCACCAGTGCCT
Reverse Sequence:	GGCGGAAGTAAGTGTGACCAC
ACCN:	<u>NM_001014979</u> , <u>NM_001014979.1</u> , <u>NM_001014979.2</u> , <u>BC092491</u> , <u>BC092491.1</u> , <u>BC042548</u> , <u>BC063391</u> , <u>BC073881</u> , <u>BC128128</u> , <u>BC128129</u> , <u>BI829598</u>
UniProt ID:	<u>A1A4V9</u>
Synonyms:	C16orf93
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)