

## Product datasheet for **HP202524**

### SLC27A6 Human qPCR Primer Pair (NM\_001017372)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	28965
Forward Sequence:	AGCCACCATGTTGTCTCACTCC
Reverse Sequence:	ACCAAAAGCCCACAGGACAGCA
ACCN:	<u>NM_001017372, NM_001017372.1, NM_001017372.2, BC041945, BX537383, NM_001017372.3</u>
UniProt ID:	<u>Q9Y2P4</u>
Synonyms:	ACSVL2; FACVL2; FATP6; VLCS-H1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)