

Product datasheet for **HP202576**

MARVELD3 Human qPCR Primer Pair (NM_001017967)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	91862
Forward Sequence:	GAAGGACTCCTGGAATGCCACA
Reverse Sequence:	GGAAAAGCTGGCACTGAATCCG
ACCN:	<u>NM_001017967</u> , <u>NM_001017967.1</u> , <u>NM_001017967.2</u> , <u>NM_001017967.3</u> , <u>BC013376</u> , <u>BC013376.2</u> , <u>BC005052</u> , <u>BC011893</u> , <u>BM742281</u> , <u>BP320419</u> , <u>BX109314</u> , <u>BX348391</u> , <u>NM_001017967.4</u>
UniProt ID:	<u>Q96A59</u>
Synonyms:	MARVD3; MRVLDC3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)