

Product datasheet for **HP202596**

RPS15 Human qPCR Primer Pair (NM_001018)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6209
Forward Sequence:	CAGCATGGTGGGCGTCTACAAC
Reverse Sequence:	GCTTGTAGGTGATGGAGAACTCG
ACCN:	<u>NM_001018</u> , <u>NM_001018.1</u> , <u>NM_001018.2</u> , <u>NM_001018.3</u> , <u>NM_001018.4</u> , <u>BC105810</u> , <u>BC105810.1</u> , <u>BC064908</u> , <u>BC080543</u> , <u>BC141832</u> , <u>BE798984</u> , <u>BG253776</u> , <u>BG705752</u> , <u>BM558707</u> , <u>BU602036</u> , <u>NM_001018.5</u>
UniProt ID:	<u>P62841</u>
Synonyms:	RIG; S15
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)