

## Product datasheet for **HP202688**

### LRR66 Human qPCR Primer Pair (NM\_001024611)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	339977
Forward Sequence:	CTGAATCCTGGAGGAAAAAGTGG
Reverse Sequence:	TGAATGGGAGGAAGCGGGTTT
ACCN:	<u><a href="#">NM_001024611</a>, <a href="#">NM_001024611.1</a>, <a href="#">NM_001024611.2</a>, <a href="#">BC040414</a>, <a href="#">BU620718</a>, <a href="#">NM_001024611.3</a></u>
UniProt ID:	<u><a href="#">Q68CR7</a></u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)