

Product datasheet for **HP202744**

ABCF1 Human qPCR Primer Pair (NM_001025091)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23
Forward Sequence:	CTCCTCAAGCACATTGCCAAC
Reverse Sequence:	GCTTCAATCGCTTGGTGTGTCAGC
ACCN:	NM_001025091 , NM_001025091.1 , BC112923 , BC016772 , BC034488 , BI547952 , BM556481 , BU620927 , NM_001025091.2
UniProt ID:	Q8NE71
Synonyms:	ABC27; ABC50
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)