

Product datasheet for **HP202811**

RPS26 Human qPCR Primer Pair (NM_001029)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6231
Forward Sequence:	GGACATTTCTGAAGCGAGCGTC
Reverse Sequence:	CGATTCCTGACTACTTTGCTGTG
ACCN:	NM_001029 , NM_001029.1 , NM_001029.2 , NM_001029.3 , BC002604 , BC002604.1 , BC013215 , BC015832 , BC070220 , BC105276 , BC105798 , NM_001029.5
UniProt ID:	P62854
Synonyms:	DBA10; eS26; S26
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)