

## Product datasheet for **HP202843**

### ARMETL1 (CDNF) Human qPCR Primer Pair (NM\_001029954)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	441549
Forward Sequence:	CAAGTCACTGATAGACAGAGGAG
Reverse Sequence:	GCTGCGTCTTTTGTGGCTCCTA
ACCN:	<a href="#">NM_001029954</a> , <a href="#">NM_001029954.2</a> , <a href="#">BC133042</a> , <a href="#">BC037872</a> , <a href="#">BC133044</a> , <a href="#">BC144619</a> , <a href="#">BG213461</a> , <a href="#">BQ025443</a> , <a href="#">NM_001029954.3</a>
UniProt ID:	<a href="#">Q49AH0</a>
Synonyms:	ARMETL1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)