

## Product datasheet for **HP202845**

### **TTL13P Human qPCR Primer Pair (NM\_001029964)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	440307
<b>Forward Sequence:</b>	GCCATCAACCTGACCAACTGCA
<b>Reverse Sequence:</b>	TGACTCGTTCCAGTGAGACAGC
<b>ACCN:</b>	<u><a href="#">NM_001029964</a>, <a href="#">NM_001029964.1</a>, <a href="#">NM_001029964.2</a>, <a href="#">BC036668</a>, <a href="#">BC144453</a>, <a href="#">BC160183</a></u>
<b>Synonyms:</b>	FLJ46079; MGC33417
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)