

Product datasheet for **HP202848**

METT11D1 (METTL17) Human qPCR Primer Pair (NM_001029991)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64745
Forward Sequence:	ACTTCGTGGAGCAGTGCTACAC
Reverse Sequence:	CAAAGCCACCATCCAGTCTTGC
ACCN:	<u>NM_001029991</u> , <u>NM_001029991.1</u> , <u>BC005053</u> , <u>BG437086</u> , <u>NM_001029991.2</u>
UniProt ID:	<u>Q9H7H0</u>
Synonyms:	METT11D1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)