

Product datasheet for **HP202926**

DDX59 Human qPCR Primer Pair (NM_001031725)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	83479
Forward Sequence:	GCCGTTTCAGAAAATCACAGGGC
Reverse Sequence:	CCAAGACTCCTGTGCTCACTAC
ACCN:	<u>NM_001031725</u> , <u>NM_001031725.1</u> , <u>NM_001031725.2</u> , <u>NM_001031725.3</u> , <u>NM_001031725.4</u> , <u>NM_001031725.5</u> , <u>BC041801</u> , <u>BC041801.1</u> , <u>BC014183</u> , <u>BC020201</u> , <u>BM683465</u> , <u>BU621052</u> , <u>BU754283</u>
UniProt ID:	<u>Q5T1V6</u>
Synonyms:	OFD5; ZNHIT5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)