

## Product datasheet for **HP203012**

### Membralin (TMEM259) Human qPCR Primer Pair (NM\_001033026)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	91304
Forward Sequence:	TCACTACCGCTTCAATGGGCAG
Reverse Sequence:	CTGAAGCAGCATCTCCTGGATG
ACCN:	<a href="#">NM_001033026</a> , <a href="#">NM_001033026.1</a> , <a href="#">BC008957</a> , <a href="#">BC008957.2</a> , <a href="#">BQ895919</a> , <a href="#">NM_001033026.2</a>
UniProt ID:	<a href="#">Q4ZIN3</a>
Synonyms:	ASBABP1; C19orf6; MBRL; MEMBRALIN; R32184_3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)