

## Product datasheet for **HP203064**

### MTUS2 Human qPCR Primer Pair (NM\_001033602)

#### Product data:

|                   |  |
|-------------------|--|
| Product Type:     | qPCR Primer Pairs  |
| Gene ID:          | 23281  |
| Forward Sequence: | GATGGTGTCTGGAATCAAGCCC   |
| Reverse Sequence: | CTGACCGATAGAAAGTTGGAGGG  |
| ACCN:             | <u><a href="#">NM_001033602</a></u> , <u><a href="#">NM_001033602.1</a></u> , <u><a href="#">NM_001033602.2</a></u> , <u><a href="#">BC032481</a></u> , <u><a href="#">BC119654</a></u> , <u><a href="#">BC119655</a></u> ,<br><u><a href="#">BC129985</a></u> , <u><a href="#">BC150245</a></u> , <u><a href="#">BM906724</a></u> , <u><a href="#">BU788127</a></u> |
| UniProt ID:       | <u><a href="#">Q5JR59</a></u>  |
| Synonyms:         | CAZIP; ICIS; KIAA0774; TIP150  |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.                         |
| Storage:          | Store at -20°C.  |
| Stability:        | The primer mix is stable for one year from date of shipping.   |



[View online »](#)