

## Product datasheet for **HP203121**

### SRA1 Human qPCR Primer Pair (NM\_001035235)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10011
Forward Sequence:	GAGGATGTGCTGAGACCTTTGG
Reverse Sequence:	TTTCCTCCAGCCCACTGTTCTT
ACCN:	<a href="#">NM_001035235</a> , <a href="#">NM_001035235.1</a> , <a href="#">NM_001035235.2</a> , <a href="#">NM_001035235.3</a> , <a href="#">BC040043</a> , <a href="#">BC067895</a> , <a href="#">BC152788</a> , <a href="#">BC156760</a> , <a href="#">BQ052761</a> , <a href="#">BU553214</a>
UniProt ID:	<a href="#">Q9HD15</a>
Synonyms:	pp7684; SRA; SRAP; STRAA1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)