

## Product datasheet for **HP203173**

### HELZ2 Human qPCR Primer Pair (NM\_001037335)

#### Product data:

|                   |   |
|-------------------|---|
| Product Type:     | qPCR Primer Pairs   |
| Gene ID:          | 85441   |
| Forward Sequence: | GAGGTGCATCTGTGTCGTTTC   |
| Reverse Sequence: | CAGGATCTCAAACTGCCGACAG  |
| ACCN:             | <u><a href="#">NM_001037335</a></u> , <u><a href="#">NM_001037335.1</a></u> , <u><a href="#">NM_001037335.2</a></u> , <u><a href="#">BC035130</a></u> , <u><a href="#">BC036883</a></u> , <u><a href="#">BC104974</a></u> , <u><a href="#">BC112023</a></u> , <u><a href="#">BC150283</a></u> , <u><a href="#">BE062855</a></u> , <u><a href="#">BE731343</a></u> , <u><a href="#">BQ707207</a></u> |
| UniProt ID:       | <u><a href="#">Q9BYK8</a></u>   |
| Synonyms:         | PDIP-1; PRIC285   |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)  |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.  |
| Storage:          | Store at -20°C.   |
| Stability:        | The primer mix is stable for one year from date of shipping.  |



[View online »](#)