

Product datasheet for **HP203220**

SLC38A10 Human qPCR Primer Pair (NM_001037984)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	124565
Forward Sequence:	ATGTGCCTGGACGGGAAGATGG
Reverse Sequence:	ACCACCTTGTCGTGAGGAACAG
ACCN:	<u>BC014156</u> , <u>NM_001037984</u> , <u>NM_001037984.1</u> , <u>NM_001037984.2</u> , <u>BC014642</u> , <u>BC014642.1</u> , <u>BC007391</u> , <u>BC011702</u> , <u>BC020925</u> , <u>BI194553</u> , <u>BM129331</u> , <u>BQ642219</u> , <u>BU552914</u> , <u>NM_001037984.3</u>
UniProt ID:	<u>Q9HBR0</u>
Synonyms:	PP1744
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)