

Product datasheet for **HP203266**

KIAA1984 (CCDC183) Human qPCR Primer Pair (NM_001039374)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84960
Forward Sequence:	GCCAGAACATCCACCTGCTGTA
Reverse Sequence:	CTCTGAGCAGTAGTTGACCACG
ACCN:	<u>NM_001039374</u> , <u>NM_001039374.1</u> , <u>NM_001039374.2</u> , <u>NM_001039374.3</u> , <u>NM_001039374.4</u> , <u>BC064582</u> , <u>BC137317</u> , <u>BC137318</u> , <u>BF979345</u> , <u>NM_001039374.5</u>
UniProt ID:	<u>Q5T5S1</u>
Synonyms:	bA216L13.7; KIAA1984; PARF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)